

OMRON


## Automatic Blood Pressure Monitor

Model HEM-716A  
Instruction Manual




Intellisense®  
3298839-8D  
Issue Date: 2024-07-12

### VOICE GUIDE FUNCTION AVAILABLE




#### ■ Language setting



5.2

#### ■ Volume setting



5.3

### Table of Contents

- Introduction
- Important Safety Information
- 1. Know Your Monitor
- 2. Preparation
- 3. Using Your Monitor
- 4. Using Memory Function
- 5. Other Settings
- 6. Error Messages and Troubleshooting
- 7. Maintenance
- 8. Optional Medical Accessories
- 9. Other Optional Part
- 10. Specifications
- 11. Guidance and Manufacturer's Declaration

### Introduction

Thank you for purchasing the OMRON Automatic Blood Pressure Monitor. This monitor can be used by multiple users. This blood pressure monitor uses the oscillometric method of blood pressure measurement. This means this monitor detects your blood movement through your brachial artery and converts the movements into a digital reading.

#### Safety Instructions

This instruction manual provides you with important information about the OMRON Automatic Blood Pressure Monitor. To ensure the safe and proper use of this monitor, READ and UNDERSTAND all of these instructions. **If you do not understand these instructions or have any questions, contact your local OMRON representative before attempting to use this monitor. For specific information about your own blood pressure, consult with your physician.**

#### Intended Use

##### Intended Purpose

This device is a digital monitor intended for use in measuring blood pressure and pulse rate in the adult patient population.

##### Intended Patients

Adult patient population

##### Intended Users

Adult population who can understand this instruction manual.

##### Clinical Benefit

Patient's blood pressure can be measured non-invasively and simply in the home environment.

#### Type of Use/Reuse

This monitor is intended to be multiple patient multiple use.

#### Limitation

Patient's arm circumference must be 17 - 42 cm.

#### Indication

This device is used by healthy individuals, patients with hypertension, patients with health-conscious individuals, in a general household situation for measuring blood pressure and pulse rate.

#### Receiving and Inspection

Remove this monitor and other components from the packaging and inspect for damage. If this monitor or any other components is damaged, DO NOT USE and consult with your local OMRON representative.

### Important Safety Information

**Read the Important Safety Information in this instruction manual before using this monitor.**

Follow this instruction manual thoroughly for your safety. Keep for future reference. For specific information about your own blood pressure, CONSULT WITH YOUR PHYSICIAN.

#### Contraindications

- DO NOT use this monitor on an injured arm or an arm under medical treatment.
- DO NOT apply the arm cuff on your arm while on an intravenous drip or blood transfusion.
- DO NOT use this monitor on infants, toddlers, children or persons who cannot express themselves.

#### Side Effects

- Taking measurements more often than necessary may cause bruising due to blood flow interference.
- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied. NOTE: refer to "If your systolic pressure is more than 210 mmHg" in section 3 for additional information.
- Stop using this monitor and consult with your physician if you experience skin irritation or discomfort.

**Warning** Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**DO NOT adjust medication (including altering your use of any drug or treatment) based on readings from this blood pressure monitor. Take medication as prescribed by your physician. ONLY a physician is qualified to diagnose and treat high blood pressure and heart related conditions.**

**NEVER diagnose or treat yourself based on your readings. ALWAYS consult with your physician.**

**If you are experiencing any symptoms or concerns, consult with your physician.**

**DO NOT postpone/stop regular checkups or physician visits based on the results that you are getting from this monitor.**

**DO NOT use this monitor in areas containing high frequency (HF) surgical equipment, magnetic resonance imaging (MRI) equipment, computerized tomography (CT) scanners. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.**

**DO NOT use this monitor in oxygen rich environments or near flammable gas.**

**Consult with your physician before using this monitor if you have common arrhythmias such as atrial or ventricular premature beats or atrial fibrillation; arterial sclerosis; poor perfusion; diabetes; pregnancy; pre-eclampsia or renal disease. NOTE that any of these conditions in addition to patient motion, trembling, or shivering may affect the measurement reading.**

**To help avoid strangulation, keep the air tube and AC adapter cable away from infants, toddlers and children.**

**This product contains small parts that may cause a choking hazard if swallowed by infants, toddlers and children.**

#### AC Adapter (optional accessory) Handling and Usage

**DO NOT use the AC adapter if this monitor or the AC adapter cable is damaged. If this monitor or the cable is damaged, turn off the power and unplug the AC adapter immediately.**

**Plug the AC adapter into the appropriate voltage outlet. DO NOT use in a multi-outlet plug.**

**NEVER plug in or unplug the AC adapter from the electric outlet with wet hands.**

**DO NOT disassemble or attempt to repair the AC adapter.**

#### Battery Handling and Usage

**Keep the batteries out of the reach of infants, toddlers and children.**

**Precaution** Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient, or cause damage to the equipment or other property.

**Consult with your physician before using this monitor on an arm where intravascular access or therapy, or an arteriovenous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury.**

**Consult with your physician before using this monitor if you have had a mastectomy or lymph node clearance.**

**Consult with your physician before using this monitor if you have severe blood flow problems or blood disorders as cuff inflation can cause bruising.**

**ONLY inflate the arm cuff when it is applied on your upper arm.**

**Remove the arm cuff if it does not start deflating during a measurement.**

**DO NOT use this monitor for any purpose other than measuring blood pressure.**

**During measurement, make sure that no mobile device or any other electrical device that emits electromagnetic fields is within 30 cm of this monitor. This may result in incorrect operation of the monitor and/or cause an inaccurate reading.**

**DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.**

**DO NOT use in a location where there is moisture or a risk of water splashing this monitor. This may damage this monitor.**

**DO NOT use this monitor in a moving vehicle such as in a car or on an aircraft.**

**DO NOT drop or subject this monitor to strong shocks or vibrations.**

**DO NOT use this monitor in places with high or low humidity or high or low temperatures. Refer to section 10.**

**During measurement, observe the arm to ensure that the monitor is not causing prolonged impairment to blood circulation.**

**DO NOT use this monitor in high-use environments such as medical clinics or physician offices.**

**DO NOT use this monitor with other medical electrical (ME) equipment simultaneously. This may result in incorrect operation of the devices and/or cause an inaccurate reading.**

**Avoid bathing, drinking alcohol or caffeine, smoking, exercising and eating for at least 30 minutes before taking a measurement.**

**Rest for at least 5 minutes before taking a measurement.**

**Remove tight-fitting or thick clothing from your arm while taking a measurement.**

**Remain still and DO NOT talk while taking a measurement.**

**ONLY use the arm cuff on persons whose arm circumference is within the specified range of the cuff.**

**Ensure that this monitor has acclimated to room temperature before taking a measurement. Taking a measurement after an extreme temperature change could lead to an inaccurate reading. It is recommended that you wait for approximately 2 hours for the monitor to warm up or cool down when the monitor is used in an environment within the temperature specified as operating conditions after it is stored either at the maximum or at the minimum storage temperature.**

**For additional information on operating and storage/transport temperature, refer to section 10.**

**DO NOT use this monitor after the durable period has ended. Refer to section 10.**

**DO NOT crease the arm cuff or the air tube excessively.**

**DO NOT fold or kink the air tube while taking a measurement. This may cause an injury by interrupting blood flow.**

**To unplug the air plug, pull on the plastic air plug at the base of the tube, not the tube itself.**

**ONLY use the AC adapter, arm cuff, batteries and accessories specified for this monitor. Use of unsupported AC adapters, arm cuffs and batteries may damage and/or may be hazardous to this monitor.**

**ONLY use the approved arm cuff for this monitor. Use of other arm cuffs may result in incorrect readings.**

#### AC Adapter (optional accessory) Handling and Usage

**Fully insert the AC adapter into the outlet.**

**When unplugging the AC adapter from the outlet, be sure to safely pull from the AC adapter. DO NOT pull from the AC adapter cable.**

**When handling the AC adapter cable:**

**DO NOT damage it.**

**DO NOT break it.**

**DO NOT tamper with it.**

**DO NOT forcibly bend or pull it.**

**DO NOT twist it.**

**DO NOT use it if it is gathered in a bundle.**

**DO NOT pinch it.**

**DO NOT place it under heavy objects.**

**Wipe any dust off of the AC adapter.**

**Unplug the AC adapter when not in use.**

**Unplug the AC adapter before wiping this monitor.**

#### Battery Handling and Usage

**DO NOT insert batteries with their polarities incorrectly aligned.**

**ONLY use 4 "AA" alkaline or manganese batteries with this monitor. DO NOT use other types of batteries. DO NOT use new and used batteries together. DO NOT use different brands of batteries together.**

**Remove the batteries if this monitor will not be used for a long period of time.**

**If battery fluid should get in your eyes, immediately rinse with plenty of clean water. Consult with your physician immediately.**

**If battery fluid should get on your skin, wash your skin immediately with plenty of clean, lukewarm water. If irritation, injury or pain persists, consult with your physician.**


**DO NOT use batteries after their expiration date.**

**Periodically check the batteries to ensure they are in good working condition.**

2. Remain still and do not talk until the entire measurement process is complete.
- The "♥" symbol flashes at every heartbeat.
- Note**
- "Cuffwrap OK" appears if the arm cuff is wrapped around the arm correctly.
  - When "Loose-fitting cuff" is displayed, the arm cuff is not applied correctly. Press the [START/STOP] button to turn your monitor off, then apply the cuff correctly.

Your blood pressure and pulse rate are displayed. The voice reads out the rates simultaneously.

If either the systolic or the diastolic reading is high (refer to sub-section 1.3), "High" appears.



3. Press the [START/STOP] button to turn your monitor off.

- Note**
- Your monitor automatically turns off after 2 minutes.
  - Wait 2 - 3 minutes between measurements. The wait time allows the arteries to return to their pre-measurement state. You may need to increase the wait time depending on your individual physiological characteristics.

- **If your systolic pressure is more than 210 mmHg**
- If your systolic pressure is more than 210 mmHg, take a measurement manually. After the arm cuff starts to inflate, press and hold the [START/STOP] button until the monitor inflates 30 to 40 mmHg higher than your expected systolic pressure.

- Note**
- Do not inflate above 299 mmHg.

**⚠ Precaution**


- Inflating to a higher pressure than necessary may result in bruising of the arm where the cuff is applied.

#### 4. Using Memory Function

Your monitor automatically stores up to 30 readings. This monitor tracks readings for 1 user. If you are sharing with another, manually track your readings separately.


##### 4.1 Viewing the Readings Stored in Memory

1. Press the [M] button.



The Memory number appears for 2 seconds before the pulse rate is displayed. The most recent reading stored in the memory will be numbered as "1".

- Note**
- If the reading was high (refer to sub-section 1.3), "High" appears.
  - If there are no readings stored in the memory, the screen to the right is displayed.
  - If the memory is full, new readings will replace the old readings, starting with the oldest reading.



2. Press the [M] button repeatedly to scroll through the previous readings stored in the memory.

**To Manage Your Readings on Your Smartphone**

To manage your readings on your smartphone, use "OMRON connect" app. Follow the steps below and keep records of several readings over a certain period of time.


- Download and install the free "OMRON connect" app onto your smartphone.
- Open the app on your smartphone and follow the instructions.
- Enter your readings to the app.

#### 5. Other Settings


- 5.1 **Restoring Your Monitor to the Default Settings**
- To delete all the information stored in your monitor, follow the instructions below. Make sure that your monitor is turned off.

1. Press the [M] button.
- After a memory number appears, the latest reading will appear.

2. While holding the [M] button down, press and hold the [START/STOP] button for more than 4 seconds.
- The display to the right flashes.



3. While holding the [M] button down again, press and hold the [START/STOP] button for more than 4 seconds.
- The display stops flashing, then your monitor is restored to the default settings.




4. Press the [START/STOP] button to turn your monitor off.

- Note**
- Your monitor automatically turns off after 10 seconds.

- 5.2 **Setting the Voice Language**
- You can choose English or Hindi for the voice guide function. English is set by default.
- To change the language, follow the instructions below. Make sure that your monitor is turned off.

1. Press the [V] button.
- The language ("En" or "Hn") and volume level you set last time flash. The voice reads out the volume simultaneously.




- Note**
- "En" stands for "English", "Hn" stands for "Hindi".

2. Press the [V] button to switch the language.
- The voice reads out the volume level in the language you set.

- Note**
- To cancel the setting, press the [START/STOP] button before confirming the setting in step 3.

3. Press and hold the [V] button for more than 3 seconds.
- The display stops flashing and the setting is confirmed.
- The voice notifies it simultaneously.




4. Press the [START/STOP] button to turn your monitor off.

- Note**
- Your monitor automatically turns off after 10 seconds.

- 5.3 **Setting the Volume**
- You can choose the volume level for the voice guide function. "1" is set as the volume level by default. To adjust the volume, follow the instructions below. Make sure that your monitor is turned off.


1. Press the [V] button.
- The volume level and the language ("En" or "Hn") you set last time flash. The voice reads out the volume simultaneously.



2. Press the [V] button repeatedly to set the volume level.
- The volume level can set between 1 to 3. If "oFF" is selected, the voice guide function will be disabled. The voice reads out the volume level in the language you set.

- Note**
- To cancel the setting, press the [START/STOP] button before confirming the setting in step 3.

3. Press and hold the [V] button for more than 3 seconds.
- The display stops flashing and the setting is confirmed.
- The voice notifies it simultaneously.



4. Press the [START/STOP] button to turn your monitor off.

- Note**
- Your monitor automatically turns off after 10 seconds.

#### 6. Error Messages and Troubleshooting

If any of the below problems occur during measurement, check to make sure that no other electrical device is within 30 cm of the monitor. If the problem persists, refer to the table below.

Display/Problem	Possible Cause	Solution
E1 appears or the arm cuff does not inflate.	The [START/STOP] button was pressed while the arm cuff is not applied.	Press the [START/STOP] button again to turn the monitor off.
	Air plug is not completely plugged into the monitor.	Insert the air plug securely.
	The arm cuff is not applied correctly.	Apply the arm cuff correctly, then take another measurement. Refer to sub-section 2.3.
E2 appears or a measurement cannot be completed after the arm cuff inflates.	Air is leaking from the arm cuff.	Replace the arm cuff to the new one. Refer to section 8.
	Moving or talking during a measurement causes the arm cuff to not inflate sufficiently.	Remain still and do not talk during a measurement. If "E2" appears repeatedly, inflate the arm cuff manually until the systolic pressure is 30 to 40 mmHg above your previous readings. Refer to sub-section 3.1.
E3 appears	The arm cuff is inflated exceeding the maximum allowable pressure.	Do not touch the arm cuff and/or bend the air tube while taking a measurement. If inflating the arm cuff manually, refer to sub-section 3.1.
E4 appears	Moving or talking during a measurement results in vibrations that disrupt the measurement.	Remain still and do not talk during a measurement.
E5 appears	Irregular heartbeat 2 times / Irregular heartbeat 3 times / Try again and Body movement / Try again and Loose-fitting cuff appears	Apply the arm cuff correctly, then take another measurement. Refer to sub-section 2.3. Remain still and sit correctly during a measurement. If "Irregular heartbeat" continues to appear, it is recommended that you consult with your physician.
does not flash during a measurement		
Er appears	The monitor has malfunctioned.	Press the [START/STOP] button again. If "Er" still appears, contact your local OMRON representative.
flashes	Batteries are low.	Replacing all 4 batteries with new ones is recommended. Refer to sub-section 2.1.
A1 and appears or the monitor is turned off unexpectedly during a measurement.	Batteries are depleted.	Immediately replace all 4 batteries with new ones. Refer to sub-section 2.1.
Batteries are depleted.		
Nothing appears on the display of the monitor.	Battery polarities are not properly aligned.	Check the battery installation for proper placement. Refer to sub-section 2.1.
No voice is heard.	The volume is set on "oFF".	Set the volume level between 1 to 3. Refer to sub-section 5.3.
You want to hear the voice in Hindi.	English is set by default. To change the voice language, refer to sub-section 5.2.	
The volume is too loud or quiet.	To change the volume, refer to sub-section 5.3.	
Readings appear too high or too low.	Blood pressure varies constantly. Many factors including stress, time of day, and/or how you apply the arm cuff, may affect your blood pressure. Review sub-sections 2.2 - 2.4 and section 3.	
Any other problems occur.	Press the [START/STOP] button to turn the monitor off, then press it again to take a measurement. If the problem continues, remove all batteries and wait for 30 seconds. Then re-install the batteries. If the problem still persists, contact your local OMRON representative.	

#### 7. Maintenance

##### 7.1 Maintenance

To protect your monitor from damage, follow the directions below:

- Changes or modifications not approved by the manufacturer will void the user warranty.

##### ⚠ Precaution

- DO NOT disassemble or attempt to repair this monitor or other components. This may cause an inaccurate reading.

##### 7.2 Storage

Keep your monitor and other components in the storage case when not in use. Store your monitor and other components in a clean, safe location.

1. Remove the arm cuff from the monitor.


The air plug is designed to fit securely into the monitor so it does not accidentally unplug from the monitor when in use. When removing the arm cuff from the monitor, pull the air plug out from the monitor in a straight motion, holding the plug and not using the tube. You will hear a "click" sound as the cuff plug is removed from the monitor.

##### ⚠ Precaution


- To unplug the air plug, pull on the plastic air plug at the base of the tube, not the tube itself.

2. Gently fold the air tube into the arm cuff.

- Note**
- Do not bend or crease the air tube excessively.



3. Place your monitor and other components in the storage case.



Do not store your monitor and other components:

- If your monitor and other components are wet.
- In locations exposed to extreme temperatures, humidity, direct sunlight, dust or corrosive vapors such as bleach.
- In locations exposed to vibrations or shocks.

##### 7.3 Wiping the Monitor

- Do not use any abrasive or volatile cleaners.
- Use a soft dry cloth or a soft cloth moistened with mild (neutral) detergent to wipe your monitor and arm cuff and then wipe them with a dry cloth.
- Do not wash or immerse your monitor and arm cuff or other components in water.
- Do not use gasoline, thinners or similar solvents to wipe your monitor and arm cuff or other components.

##### 7.4 Calibration and Service


- The accuracy of this blood pressure monitor has been carefully tested and is designed for a long service life.
- It is generally recommended to have the unit inspected every two years to ensure correct functioning and accuracy. Please consult your authorised OMRON dealer or the OMRON Customer Service at the address given on the packaging or attached literature.

#### 8. Optional Medical Accessories



- Note**
- Do not throw the air plug away. The air plug can be applicable to the optional cuff.


#### 9. Other Optional Part



##### Using the AC Adapter (optional accessory)

- Note**
- Make sure not to place your monitor in a location where it is difficult to plug and unplug the AC adapter.

1. Insert the AC adapter plug into the AC adapter jack on the back side of your monitor.



2. Plug the AC adapter into an electrical outlet.
- To unplug the AC adapter, unplug the AC adapter from the electrical outlet and then remove the AC adapter plug from the monitor.

#### ⚠ Precaution

- ONLY use the AC adapter, arm cuff, batteries and accessories specified for this monitor. Use of unsupported AC adapters, arm cuffs and batteries may damage and/or may be hazardous to this monitor.

#### 10. Specifications

Product description	Automatic Blood Pressure Monitor
Model	HEM-716A
Display	LCD digital display
Cuff pressure range	0 to 299 mmHg
Blood pressure measurement range	SYS: 60 to 260 mmHg DIA: 40 to 215 mmHg
Pulse measurement range	40 to 180 beats / min.
Accuracy	Pressure: ±3 mmHg Pulse: ±5 % of display reading
Inflation	Automatic by electric pump
Deflation	Automatic pressure release valve
Measurement method	Oscillometric method
Operating mode	Continuous operation
IP classification	Monitor: IP20 Optional AC adapter: IP22
Rating	DC6 V 4 W
Power source	4 "AA" batteries 1.5 V or optional AC adapter (INPUT AC 100 - 240 V 50 - 60 Hz 0.12 - 0.065 A)
Battery life	Approximately 500 measurements (using new alkaline batteries and included arm cuff. Depending on the type of battery and arm cuff size.) Monitor: 5 years or the time when reaching 10,000 times of use. Arm cuff: 5 years or the time when reaching 10,000 times of use. Optional AC adapter: 5 years
Durable period (Service life)	+10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Operating conditions	+10 to +40 °C / 15 to 90 % RH (non-condensing) / 800 to 1060 hPa
Storage / Transport conditions	-20 to +60 °C / 10 to 90 % RH (non-condensing)
Weight	Monitor: approximately 340 g (not including batteries) Arm cuff: approximately 170 g Monitor: Approximately 109 mm (w) × 83 mm (h) × 134 mm (l) Arm cuff: Approximately 145 mm × 594 mm (air tube: 750 mm)
Dimensions	Monitor: Approximately 109 mm (w) × 83 mm (h) × 134 mm (l) Arm cuff: Approximately 145 mm × 594 mm (air tube: 750 mm)
Cuff circumference applicable to the monitor	Total arm circumference range: 17 to 42 cm Included arm cuff: 22 to 42 cm Optional arm cuff: 17 to 22 cm
Memory	Stores up to 30 readings
Contents	Monitor, arm cuff (HEM-RML31 Type B), 4 "AA" batteries, storage case, instruction manual, warranty card
Protection against electric shock	Internally powered ME equipment (When using only batteries) Class II ME equipment (Optional AC adapter)
Applied part	Type BF (arm cuff)

##### Note

- These specifications are subject to change without notice.
- This monitor has not been validated for use on pregnant patients.
- IP classification is degrees of protection provided by enclosures in accordance with IEC 60529. This monitor and optional AC adapter are protected against solid foreign objects of 12.5 mm diameter and greater such as a finger. The optional AC adapter is protected against oblique falling water drops which may cause issues during a normal operation.
- Operation mode classification complies with IEC 60601-1.

#### 11. Guidance and Manufacturer's Declaration

- This OMRON product is produced under the strict quality system of OMRON HEALTHCARE Co., Ltd., Japan. The Core component for OMRON blood pressure monitors, which is the Pressure Sensor, is produced in Japan.
- Please report to the manufacturer and the regulatory authority in which you are established about any serious incident that has occurred in relation to this device.

Symbols description	
	Applied part - Type BF Degree of protection against electric shock (leakage current)
	CE Marking
	Class II equipment. Protection against electric shock
IP XX	Ingress protection degree provided by IEC 60529
	Serial number
	LOT number
	Medical Device
	Temperature limitation

Symbols description	
	Humidity limitation
	Atmospheric pressure limitation
	Indication of connector polarity
	Identifier of cuffs compatible for the device
	OMRON's trademarked technology for blood pressure measurement
	Arm circumference
	Cuff positioning indicator for the left arm
	Artery Mark
	Range pointer and brachial artery alignment position
	Range indicator of arm circumferences to help selection of the correct cuff size
	Need for the user to consult this instruction manual.
	Need for the user to follow this instruction manual thoroughly for your safety. (Background: blue)
	Direct current
	Alternating current
	Date of manufacture
	BIS mark
	Efficiency level of power supply
	Not made with natural rubber latex
	Recycling mark
	Catalogue number

#### Important information regarding Electromagnetic Compatibility (EMC)


This device conforms to the EN 60601-1-2:2015+A1:2021 Electromagnetic Compatibility (EMC) standard.

Further documentation in accordance with this EMC standard is available at <https://healthcare.omron.com/electro-magnetic-compatibility>

#### Correct Disposal of This Product (Waste Electrical & Electronic Equipment)

This marking shown on the product or its literature, indicates that it should not be disposed of, with other household wastes at the end of its working life.

To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate this product from other types of wastes and recycle it responsibly to promote the sustainable reuse of material resources.



Household users should contact either the retailer where they purchased this product, or their local government office, for details of where and how they can return this item for environmentally safe recycling.

Business users should contact their supplier and check the terms and conditions of the purchase contract. This product should not be mixed with other commercial waste for disposal.

Apple and the Apple logo are trademarks of Apple Inc., registered in the U.S. and other countries and regions. App Store is a service mark of Apple Inc.

Google Play and the Google Play logo are trademarks of Google LLC.

Other trademarks and trade names are those of their respective owners.

<https://www.omronhealthcare-ap.com/>

<b>Manufacturer</b>	OMRON HEALTHCARE Co., Ltd. 53, Kunotsubo, Terado-cho, Muko, KYOTO, 617-0002 JAPAN
<b>EU-representative</b>	OMRON HEALTHCARE EUROPE B.V. Wegalaan 73, 2132 JD Hoofddorp, THE NETHERLANDS <a href="http://www.omron-healthcare.com">www.omron-healthcare.com</a>
<b>Importer in EU</b>	
<b>Marketed and Distributed by</b>	OMRON HEALTHCARE INDIA Pvt. Ltd. 6th Floor, B-Block, Sewa Tower, Plot No-19, Sector-18, Udyog Vihar, Maruti Industrial Complex, Gurugram, Haryana (India) – 122008 <a href="http://www.omronhealthcare-ap.com/in">www.omronhealthcare-ap.com/in</a>
<b>Production Facility</b>	OMRON HEALTHCARE MANUFACTURING VIETNAM CO., LTD. No.28 VSIP II, Street 2, Vietnam-Singapore Industrial Park II, Binh Duong Industry-Service-Urban Complex, Hoa Phu Ward, Thu Dau Mot City, Binh Duong Province, Vietnam